

**Safety Data Sheet**  
acc. to OSHA HCS

Printing date 11/08/2022

Reviewed on 11/02/2022

**1 Identification**

- **Product identifier**
- **Trade name:** Opal™ Bond™ Flow & Opal™ Bond™ Flow Blue
- **Article number:** SDS 75-001.12, 71161, 71050
- **Application of the substance / the mixture** Professional Orthodontic Adhesive
- **Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**  
Ultradent Products Inc.  
505 W. Ultradent Drive (10200 S)  
South Jordan, UT 84095-3942  
USA  
onlineordersupport@ultradent.com
- **Information department:** Customer Service
- **Emergency telephone number:**  
CHEMTREC (NORTH AMERICA) : (800) 424-9300  
(INTERNATIONAL) : +(703) 527-3887

**2 Hazard(s) identification**

- **Classification of the substance or mixture**



- Skin Irritation 2     H315 Causes skin irritation.
- Eye Irritation 2A     H319 Causes serious eye irritation.
- Sensitization - Skin 1     H317 May cause an allergic skin reaction.

- **Label elements**
- **GHS label elements** Void
- **Hazard pictograms** GHS07
- **Signal word** Warning
- **Health Hazard-determining components of labeling:**  
Triethylene Glycol Dimethacrylate  
Organophosphine Oxide
- **Hazard statements**  
H315 Causes skin irritation.  
H319 Causes serious eye irritation.  
H317 May cause an allergic skin reaction.
- **Precautionary statements**  
P261     Avoid breathing dust/fume/gas/mist/vapors/spray  
P264     Wash thoroughly after handling.  
P272     Contaminated work clothing must not be allowed out of the workplace.  
P280     Wear protective gloves / eye protection / face protection.  
P302+P352     If on skin: Wash with plenty of water.  
P305+P351+P338     If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.  
P362+P364     Take off contaminated clothing and wash it before reuse.

(Contd. on page 2)

# Safety Data Sheet

acc. to OSHA HCS

Printing date 11/08/2022

Reviewed on 11/02/2022

**Trade name: Opal™ Bond™ Flow & Opal™ Bond™ Flow Blue**

(Contd. of page 1)

- P333+P313 *If skin irritation or rash occurs: Get medical advice/attention.*  
 P321 *Specific treatment (see on this label).*  
 P337+P313 *If eye irritation persists: Get medical advice/attention.*  
 P363 *Wash contaminated clothing before reuse.*  
 P501 *Dispose of contents/container in accordance with local/regional/national/international regulations.*

· **Classification system:**

· **NFPA ratings (scale 0 - 4)**



· **HMIS-ratings (scale 0 - 4)**

HEALTH	1	Health = 1
FIRE	0	Fire = 0
REACTIVITY	0	Reactivity = 0

### 3 Composition/information on ingredients

- **Chemical characterization: Mixtures**  
 · **Description:** Mixture of the substances listed below with nonhazardous additions.

· **Dangerous components:**

109-16-0	Triethylene Glycol Dimethacrylate	>1-<20%
	Trade Secret	>0.1-<10%
13463-67-7	Titanium Dioxide	>0.1-<10%
162881-26-7	Organophosphine Oxide	>0.1-<10%

· **Additional information:**

The specific chemical identity of composition is being withheld as a trade secret. The specific chemical identity is made available to health professionals, employees, and designated representatives in accordance with the applicable provisions of paragraph §1910.1200.

### 4 First-aid measures

- **Description of first aid measures**  
 · **General information:** Immediately remove any clothing soiled by the product.  
 · **After inhalation:**  
*This product is a thick paste, therefore inhalation is extremely unlikely.*  
*Supply fresh air and to be sure call for a doctor.*  
*In case of unconsciousness place patient stably in side position for transportation.*  
 · **After skin contact:** Immediately wash with water and soap and rinse thoroughly.  
 · **After eye contact:**  
*Rinse opened eye for several minutes under running water. If symptoms persist, consult a doctor.*  
 · **After swallowing:** If swallowed in large quantities seek medical attention.  
 · **Information for doctor:**  
 · **Most important symptoms and effects, both acute and delayed** No further relevant information available.

(Contd. on page 3)

US

# Safety Data Sheet

## acc. to OSHA HCS

Printing date 11/08/2022

Reviewed on 11/02/2022

**Trade name: Opal™ Bond™ Flow & Opal™ Bond™ Flow Blue**

(Contd. of page 2)

- **Indication of any immediate medical attention and special treatment needed**  
No further relevant information available.

### 5 Fire-fighting measures

- **Extinguishing media**
- **Suitable extinguishing agents:**  
Foam, dry chemical, carbon dioxide  
Use fire fighting measures that suit the environment.
- **Special hazards arising from the substance or mixture** No further relevant information available.
- **Advice for firefighters**
- **Protective equipment:** Wear fully protective suit.

### 6 Accidental release measures

- **Personal precautions, protective equipment and emergency procedures** Wear protective clothing.
- **Environmental precautions:** Do not allow to enter sewers/ surface or ground water.
- **Methods and material for containment and cleaning up:**  
Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).  
Dispose contaminated material as waste according to item 13.  
Ensure adequate ventilation.
- **Reference to other sections**  
See Section 7 for information on safe handling.  
See Section 8 for information on personal protection equipment.  
See Section 13 for disposal information.

### 7 Handling and storage

- **Handling:**
- **Precautions for safe handling**  
Avoid release to the environment  
Ensure good ventilation/exhaustion at the workplace.  
Prevent formation of aerosols.
- **Information about protection against explosions and fires:** No special measures required.
- **Conditions for safe storage, including any incompatibilities**
- **Storage:**
- **Requirements to be met by storerooms and receptacles:** No special requirements.
- **Information about storage in one common storage facility:** Not required.
- **Further information about storage conditions:**  
See product labelling.  
Keep receptacle tightly sealed.
- **Specific end use(s)** Professional Orthodontic Adhesive

### 8 Exposure controls/personal protection

- **Additional information about design of technical systems:** No further data; see item 7.
- **Control parameters**
- **Components with limit values that require monitoring at the workplace:**  
The following constituent is the only constituent of the product which has a PEL, TLV or other recommended exposure limit.

(Contd. on page 4)

# Safety Data Sheet

acc. to OSHA HCS

Printing date 11/08/2022

Reviewed on 11/02/2022

**Trade name: Opal™ Bond™ Flow & Opal™ Bond™ Flow Blue**

(Contd. of page 3)

*At this time, the other constituents have no known exposure limits.*

### 13463-67-7 Titanium Dioxide

ACGIH TLV	Short-term value: 10* 5 mg/m <sup>3</sup>
PEL	Long-term value: 15* mg/m <sup>3</sup> *total dust
REL	See Pocket Guide App. A
TLV	Long-term value: 0.2* 2.5** mg/m <sup>3</sup> resp. fraction, *nanoscale, **finescale, A3
TWA	Short-term value: 15* 5 mg/m <sup>3</sup>

· **Additional information:** The lists that were valid during the creation were used as basis.

· **Exposure controls**

· **Personal protective equipment:**

· **General protective and hygienic measures:**

Keep away from foodstuffs, beverages and feed.

Immediately remove all soiled and contaminated clothing.

Wash hands before breaks and at the end of work.

Avoid contact with the eyes and skin.

· **Breathing equipment:**

In case of brief exposure or low pollution use respiratory filter device. In case of intensive or longer exposure use respiratory protective device that is independent of circulating air.

· **Protection of hands:**



Protective gloves

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.

Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.

Selection of the glove material is based on consideration of the penetration times, rates of diffusion and the degradation

· **Material of gloves**

The selection of suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

· **Penetration time of glove material**

The exact breakthrough time has to be found out by the manufacturer of the protective gloves and has to be observed.

· **Eye protection:**



Tightly sealed goggles

· **Body protection:** Protective work clothing

US

(Contd. on page 5)

# Safety Data Sheet

acc. to OSHA HCS

Printing date 11/08/2022

Reviewed on 11/02/2022

Trade name: *Opal™ Bond™ Flow & Opal™ Bond™ Flow Blue*

(Contd. of page 4)

## 9 Physical and chemical properties

· <b>Information on basic physical and chemical properties</b>	
· <b>General Information</b>	
· <b>Appearance:</b>	
· <b>Form:</b>	Paste
· <b>Color:</b>	According to product specification
· <b>Odor:</b> Acrylic	
· <b>Odor threshold:</b> Not determined.	
· <b>pH-value:</b> Not applicable (non-aqueous)	
· <b>Change in condition</b>	
· <b>Melting point/Melting range:</b>	Undetermined.
· <b>Boiling point/Boiling range:</b>	Undetermined
· <b>Flash point:</b> Not applicable	
· <b>Flammability (solid, gaseous):</b> Not applicable.	
· <b>Decomposition temperature:</b> Not determined.	
· <b>Auto igniting:</b> Product is not selfigniting.	
· <b>Danger of explosion:</b> Product does not present an explosion hazard.	
· <b>Explosion limits:</b>	
· <b>Lower:</b>	Not determined.
· <b>Upper:</b>	Not determined.
· <b>Vapor pressure:</b> Not determined.	
· <b>Density:</b> Not determined	
· <b>Relative density</b> Not determined	
· <b>Vapor density</b> Not determined.	
· <b>Evaporation rate</b> Not determined.	
· <b>Solubility in / Miscibility with</b>	
· <b>Water:</b>	Not miscible or difficult to mix.
· <b>Partition coefficient (n-octanol/water):</b> Not determined.	
· <b>Viscosity:</b>	
· <b>Dynamic:</b>	Not determined.
· <b>Kinematic:</b>	Not determined
· <b>Other information</b> No further relevant information available.	

## 10 Stability and reactivity

- **Reactivity Stable**
- **Chemical stability**
- **Thermal decomposition / conditions to be avoided:** No decomposition if used according to specifications.
- **Possibility of hazardous reactions** No dangerous reactions known.
- **Conditions to avoid**
  - Light
  - Sparks
  - Ignition sources
  - Heat
  - Flames

(Contd. on page 6)

# Safety Data Sheet

acc. to OSHA HCS

Printing date 11/08/2022

Reviewed on 11/02/2022

**Trade name: Opal™ Bond™ Flow & Opal™ Bond™ Flow Blue**

(Contd. of page 5)

- **Incompatible materials:** No further relevant information available.
- **Hazardous decomposition products:** No dangerous decomposition products known.

## 11 Toxicological information

- **Information on toxicological effects**
- **Acute toxicity:**

- **LD/LC50 values that are relevant for classification:**

### 109-16-0 Triethylene Glycol Dimethacrylate

Oral	LD50	>5,000 mg/kg (rat)
	LC50 Fish	16.4 mg/l (Fish) (Toxicity to fish)
Dermal	LD50	>2,000 mg/kg (mouse)

### 13463-67-7 Titanium Dioxide

Oral	LD50	>5,000 mg/kg (rat)
Dermal	LD50	>5,000 mg/kg (rabbit)

### 162881-26-7 Organophosphine Oxide

Oral	LD50	>2,000 mg/kg (rat)
	LC50 Fish	>0.09 mg/l (Fish) (Toxicity to fish)
Dermal	LD50	>2,000 mg/kg (rat)

- **Primary irritant effect:**
- **on the skin:** Irritant to skin and mucous membranes.
- **on the eye:** Irritating effect.
- **Sensitization:** Sensitization possible through skin contact.

- **Additional toxicological information:**

The product shows the following dangers according to internally approved calculation methods for preparations:  
Irritant

- **Carcinogenic categories**

- **IARC (International Agency for Research on Cancer)**

	Co-Cr-AL Spinel Blue Green	2B
14808-60-7	Silica Glass	I

- **NTP (National Toxicology Program)**

14808-60-7	Silica Glass	K
------------	--------------	---

- **OSHA-Ca (Occupational Safety & Health Administration)**

None of the ingredients is listed.

## 12 Ecological information

- **Toxicity**

- **Aquatic toxicity:**

### 109-16-0 Triethylene Glycol Dimethacrylate

EC50	>100 mg/kg (Algae)
Biodegradability	28 days (Aerobic) (Biodegradability testing)
Aqua toxicity	32 mg/l (daphnia) (No Observed Effect Concentration)

(Contd. on page 7)

# Safety Data Sheet

acc. to OSHA HCS

Printing date 11/08/2022

Reviewed on 11/02/2022

**Trade name: Opal™ Bond™ Flow & Opal™ Bond™ Flow Blue**

(Contd. of page 6)

### 13463-67-7 Titanium Dioxide

EC50	>100 mg/kg (Algae)
	>1,000 mg/kg (Fish)

### 162881-26-7 Organophosphine Oxide

EC50 (static)	>1.175 mg/kg (daphnia) (Toxicity to aquatic invertebrates)
Aqua toxicity	≥0.008 mg/l (daphnia) (Daphnia Magna Reproduction Test)
Toxicity to Aquatic Plants (static)	>0.26 mg/l (Plant) (Toxicity to algae)

- **Persistence and degradability** No further relevant information available.
- **Behavior in environmental systems:**
- **Bioaccumulative potential** No further relevant information available.
- **Mobility in soil** No further relevant information available.
- **Additional ecological information:**
- **General notes:**  
Water hazard class 1 (Self-assessment): slightly hazardous for water  
Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.
- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.
- **Other adverse effects** No further relevant information available.

## 13 Disposal considerations

- **Waste treatment methods**
- **Recommendation:**  
Dispose of contents/container in accordance with international, federal, state, and local regulations.
- **Uncleaned packagings:**
- **Recommendation:** Disposal must be made according to official regulations.

## 14 Transport information

- |  |                 |
|--|-----------------|
| · <b>UN-Number</b>   |                 |
| · <b>DOT, ADN, IMDG, IATA</b>  | Not Regulated   |
| · <b>UN proper shipping name</b>   |                 |
| · <b>DOT, ADN, IMDG, IATA</b>  | Not Regulated   |
| · <b>Transport hazard class(es)</b>  |                 |
| · <b>DOT, ADN, IMDG, IATA</b>  |                 |
| · <b>Class</b>   | Not Regulated   |
| · <b>Packing group</b>   |                 |
| · <b>DOT, IMDG, IATA</b>   | Not Regulated   |
| · <b>Environmental hazards:</b>  | Not Applicable. |
| · <b>Special precautions for user</b>  | Not Applicable  |
| · <b>Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code</b> | Not Applicable. |

(Contd. on page 8)

US

# Safety Data Sheet

acc. to OSHA HCS

Printing date 11/08/2022

Reviewed on 11/02/2022

**Trade name: Opal™ Bond™ Flow & Opal™ Bond™ Flow Blue**

(Contd. of page 7)

· **UN "Model Regulation":** Not Regulated

## 15 Regulatory information

· **Safety, health and environmental regulations/legislation specific for the substance or mixture**  
No further relevant information available.

· **Sara**

· **Section 355 (extremely hazardous substances):**

None of the ingredients is listed.

· **Section 313 (Specific toxic chemical listings):**

Trade Secret

Co-Cr-AL Spinel Blue Green

· **TSCA (Toxic Substances Control Act):**

109-16-0	Triethylene Glycol Dimethacrylate	ACTIVE
13463-67-7	Titanium Dioxide	ACTIVE
162881-26-7	Organophosphine Oxide	ACTIVE

· **Hazardous Air Pollutants**

Co-Cr-AL Spinel Blue Green

· **Proposition 65**

· **Chemicals known to cause cancer:**

14808-60-7 Silica Glass

· **Chemicals known to cause reproductive toxicity for females:**

None of the ingredients is listed.

· **Chemicals known to cause reproductive toxicity for males:**

None of the ingredients is listed.

· **Chemicals known to cause developmental toxicity:**

None of the ingredients is listed.

· **Carcinogenic categories**

· **EPA (Environmental Protection Agency)**

None of the ingredients is listed.

· **ACGIH Carcinogenicity (American Conference of Governmental Industrial Hygienists)**

	Trade Secret	A4
14808-60-7	Silica Glass	A2

· **NIOSH-Ca (National Institute for Occupational Safety and Health)**

14808-60-7 Silica Glass

· **Chemical safety assessment:**

Device is biocompatible when used as directed by dental professionals per ISO 10993-1

## 16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

· **Department issuing SDS:** Environmental, Health, and Safety

(Contd. on page 9)

# Safety Data Sheet

acc. to OSHA HCS

Printing date 11/08/2022

Reviewed on 11/02/2022

**Trade name: Opal™ Bond™ Flow & Opal™ Bond™ Flow Blue**

(Contd. of page 8)

- **Contact:** Customer Service
- **Date of preparation / last revision** 11/08/2022 / -
- **Abbreviations and acronyms:**
  - IMDG: International Maritime Code for Dangerous Goods
  - DOT: US Department of Transportation
  - IATA: International Air Transport Association
  - EINECS: European Inventory of Existing Commercial Chemical Substances
  - ELINCS: European List of Notified Chemical Substances
  - CAS: Chemical Abstracts Service (division of the American Chemical Society)
  - NFPA: National Fire Protection Association (USA)
  - HMS: Hazardous Materials Identification System (USA)
  - LC50: Lethal concentration, 50 percent
  - LD50: Lethal dose, 50 percent
  - PBT: Persistent, Bioaccumulative and Toxic
  - vPvB: very Persistent and very Bioaccumulative
  - NIOSH: National Institute for Occupational Safety
  - OSHA: Occupational Safety & Health
  - TLV: Threshold Limit Value
  - PEL: Permissible Exposure Limit
  - REL: Recommended Exposure Limit
  - Skin Irritation 2: Skin corrosion/irritation – Category 2
  - Eye Irritation 2A: Serious eye damage/eye irritation – Category 2A
  - Sensitization - Skin 1: Skin sensitisation – Category 1
- **\* Data compared to the previous version altered.**

US